

MAGNUM L

MAGNETOTHERAPY | 8 PROGRAMS

PHYSIOTHERAPY PRODUCT SUITABLE TO RECOVER THE FUNCTIONALITY COMPROMISED BY TRAUMA.

TECHNICAL FEATURES

Display

Visible area size 2,6"

Channels

1

Emitted intensity

120 average Gauss

Adjustable frequency

5-200 Hz

Duty Cycle

5-50%

Power supply

Mains

Size

100x160x35 mm


Weight

260 gr

EQUIPMENT

Code **G3947**

- 1 Bag
- 1 MAGNUM L magnetotherapy unit
- 1 "Flexible" solenoid, 300x100 mm
- 1 Elastic band with velcro 800x80 mm
- 1 Power supply unit
- 1 Operating manual



CE
0476



- LAST 10
- RENTAL OPTION
- MAINS SUPPLY

GLOBUS
ITALIAN EXCELLENCE

MAGNUM L

MAGNETOTHERAPY | 8 PROGRAMS

PROGRAMS LIST

HEALTH and REHAB

Small fractures

Big fractures

Cervical Osteoarthritis

Back Pain

Osteoporosis

Arthrosis

Arthritis

Epicondylitis