

DIACARE 7000

TECARTHERAPY

DIACARE 7000 IS THE MOST POWERFUL, VERSATILE AND COMPLETE DEVICE IN THE GLOBUS TECAR THERAPY RANGE. IT USES BOTH THE RESISTIVE AND THE CAPACITIVE MODE FOR THE TREATMENT OF ACUTE AND CHRONIC PATHOLOGIES OF THE MUSCULAR AND SKELETAL SYSTEM. THE DEVICE CAN BE USED ON THE PATIENT BOTH IN STATIC MANUAL MODE AND COMBINED WITH PASSIVE MOVEMENTS, OR, THANKS TO THE HANDS-FREE KIT, IT CAN BE USED DURING ACTIVE MOVEMENTS OF THE SUBJECT. DIACARE 7000 CAN ALSO BE USED WITH SPECIFIC ERGONOMIC TOOLS DESIGNED FOR THE MOBILIZATION AND THE ADVANCED TREATMENT OF ALL SOFT TISSUE DISEASES (FASCIA TOOLS).

TECHNICAL FEATURES

Display

Size 7" color touch screen display

Nominal power (NP)

350W \pm 10%

Emission frequency

300kHz, 400kHz, 448kHz, 470kHz, 700kHz, 1000kHz, 1200kHz, \pm 10%

Treatment modality

Capacitive and resistive

Capacitive electrodes

3 size \varnothing (mm 30, 50, 70)

Resistive electrodes

3 size \varnothing (mm 30, 50, 70)

Mains supply

- 230V AC, 50-60Hz, 3.7A
- 115V AC, 50-60Hz, 7.4A

Size

515x340x170 mm

Weight

11,5 Kg

EQUIPMENT

Cod. G5210

- 1 DIACARE 7000 tecartherapy unit
- 1 Counter electrode (plate)
- 1 Monopolar resistive handpiece with 3 heads \varnothing (mm 30, 50, 70)
- 1 Monopolar capacitive handpiece with 3 heads \varnothing (mm 30, 50, 70)
- 1 Bipolar headpiece \varnothing (mm 50)
- 1 Kit for automatic treatment
- 1 Cream can ml 1000
- 1 Power supply unit
- 1 Operating manual



DIACARE 7000

TECARTHERAPY

PRE-SET PROGRAMS TO TREAT THE MOST COMMON PATHOLOGIES

BONE AND ARTICULAR PATHOLOGIES

PATHOLOGIES OF TENDONS AND LIGAMENTS

MUSCULAR AND SOFT TISSUE PATHOLOGIES

VASCULAR AND LYMPHATIC PATHOLOGIES

PATHOLOGIES OF THE NERVOUS SYSTEM



ATLAS



PULSED TREATMENTS



ATHERMAL TREATMENTS



MULTI-FREQUENCY



AUTOMATIC FREQUENCY SCANNING



MYOFASCIALS



HANDS-FREE MODE



TRANSDERMAL VEHICULATION



RESISTIVE



CAPACITIVE



BIPOLAR



COLOR DISPLAY



TOUCH SCREEN DISPLAY



MAINS SUPPLY