

# MAGNUM 2500

## MAGNETOTHERAPY | 52 PROGRAMS

PORTABLE MAGNETOTHERAPY, IDEAL FOR HOME TREATMENTS. MAGNUM 2500 IS RECOMMENDABLE FOR A BETTER RECOVERY FROM TRAUMAS AND FRACTURES, BUT ALSO TO GIVE RELIEF AND HELP TO ARTICULAR REGENERATION IN CHRONIC PATHOLOGIES SUCH AS ARTHROSIS AND ARTHRITIS.

### TECHNICAL SPECIFICATIONS

#### Display

Size of the visible area 2.6"

#### Channels

2

#### Intensity can be emitted with flexible diffuser

- 160 Gauss on average for channel
- 320 Gauss on average on the whole (for 2 channels)

#### Frequency

5-200 Hz

#### Duty Cycle

5-50%

#### Free programming

Possibility of adding up to 5 custom programs

#### 1+1 Function

Possibility of running 2 different programs at the same time, using 1 channel each

#### Power supply

- Mains
- Rechargeable batteries (not included)

#### Size

100x160x35 mm


#### Weight

280 gr

### STANDARD EQUIPMENT

Code **G5438**

- 1 Bag
- 1 MAGNUM 2500 magnetotherapy device
- 2 "Flexibles" solenoids, 300x100 mm
- 2 Elastic bands with Velcro 800x80 mm
- 1 Magnum Test
- 1 Power supply unit
- 1 Operating manual



### PURCHASE OPTION

Code **G6191** MAGNUM 2500 with 2 solenoids "Pocket Pro" and 1 elastic band with Velcro



**CE**  
0476



- EASY PROGRAM
- 1+1 MODE
- LAST 10
- MULTI USER
- PROGRAMMABLE
- RENTAL OPTION
- RUN TIME / EDIT PARAMETERS
- MAINS SUPPLY

**GLOBUS**  
ITALIAN EXCELLENCE

# MAGNUM 2500

## MAGNETOTHERAPY | 52 PROGRAMS

### PROGRAM LIST

#### ARTHROSIS/ARTHRITIS

Knee arthrosis

Cervical arthrosis

Arthrotos

Arthrosis

Hand arthrosis

Hip arthrosis

Shoulder arthrosis

Column arthrosis

#### PAIN

Algodystrophy - Sudeck disease

Back pain

Joint pain

#### INFLAMMATION

Anti-inflammatory

Epicondylitis

Epitrochleytis

Carpal tunnel

Shoulder tendinitis

#### ORTHOPEDICS

Posto ACL surgery

Cartilage lesions

Knee prothesis

Hip prothesis

#### BONE PATHOLOGIES

Fractures

Scaphoid fracture

Wrist fracture

Shoulder Fracture

Collarbone fracture

Forearm fracture

Pelvic fracture

Tibial fracture

Femoral fracture

Delayed fracture consolidation

Fracture with plaster

Osteonecrosis

Osteoporosis

Bone edemas

#### TRAUMAS / EDEMAS / WOUNDS

Ankle sprain

Edemas- contusion

#### DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY

Multifrequency program

Program 10 Hz

Program 20 Hz

Program 30 Hz

Program 40 Hz

Program 50 Hz

Program 60 Hz

Program 70 Hz

Program 80 Hz

Program 90 Hz

Program 100 Hz

Program 120 Hz

Program 140 Hz

Program 160 Hz

Program 180 Hz

Program 200 Hz