

MAGNUM XL PRO

MAGNETOTHERAPY | 46 PROGRAMS

MAGNETOTHERAPY DEVICE EQUIPPED WITH A FULL RANGE OF PRE-SET PROGRAMS TO TREAT THE MOST COMMON OSTEOARTICULAR DISORDERS. THE DEVICE IS ONE OF THE MOST POWERFUL ON THE MARKET AND IT IS PROVIDED WITH TWO SOLENOIDS.

TECHNICAL FEATURES

Display

Visible area size 2.6"

Channels

2

Emitted intensity

- 200 average Gauss per channel
- 400 Gauss on average on the whole (for 2 channels)

Adjustable frequency

5-200 Hz

Duty Cycle

5-50%

Free programming

Possibility of adding up to 15 custom programs

Power supply

- Mains
- Rechargeable batteries (not included)

Size

100x160x35 mm


Weight

280 gr

EQUIPMENT

Code **G3956**

- 1 Bag
- 1 MAGNUM XL PRO magnetotherapy unit
- 2 "Flexibles" solenoids, 300x100 mm
- 2 Elastic bands with Velcro 800x80 mm
- 1 Magnum Test
- 1 Power supply unit
- 1 Operating Manual



PURCHASE OPTION

Code **G6033** MAGNUM XL PRO with 2 solenoids "Pocket Pro" and 1 elastic band with Velcro



CE
0476



- EASY PROGRAM
- LAST 10
- PROGRAMMABLE
- MULTI USER
- RENTAL OPTION
- RUN TIME / EDIT PARAMETERS
- MAINS SUPPLY

GLOBUS
ITALIAN EXCELLENCE

MAGNUM XL PRO

MAGNETOTHERAPY | 46 PROGRAMS

PROGRAMS LIST

HEALT and REHAB

Fractures

Shoulder Fractures

Collarbone Fractures

Forearm Fractures

Hip Fractures

Tibial Fractures

Wrist Fracture

Femoral neck fracture

Scaphoid Fracture

Delayed Fracture Consolidation

Rotator Cuffs

Carpal Tunnel

Back Pain

Arthrosis

Osteoporosis

Cervical Osteoarthritis

Hand Osteoarthritis

Shoulder Osteoarthritis

Spinal Osteoarthritis

Knee Osteoarthritis

Hip Osteoarthritis

Hip Prosthesis

Post ACL Surgery

Knee Clean-out

Knee Prosthesis

Epitrochleitis

Epicondylitis

Tendinitis

Edemas - Confusions

Articular pains

Ankle sprain

Arthritis

BEAUTY and FITNESS

Abdomen Drainage

Upper Limbs Drainage

Thigh, Gluteus, Hip Drain

Abdomen Cellulite

Upper Limbs Cellulite

Thigh, Gluteus and Hip cellulite

Oedematous States

Swollen Legs

Reactivating Treatment

Tissue Regeneration

Acne

Scars

Face Wrinkles

Stretch Marks