

TECARVET 4000

TECAR THERAPY

TECARVET 4000 IS A DEVICE THAT USES DIATHERMY TO OBTAIN EXTRAORDINARY THERAPEUTIC EFFECTS ON PETS SUCH AS CATS AND DOGS. THE DEVICE TRANSFERS ENERGY IN BOTH CAPACITIVE AND RESISTIVE MODE. IN THIS WAY THE EFFECT OF THE TREATMENT CAN BE DIRECTED WITH GOOD EFFECTIVENESS BOTH ON DEEP BIOLOGICAL TISSUES AND ON SUPERFICIAL ONES.

SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

TECHNICAL FEATURES

Display

Size 4.3" touch-screen

Power

160W, $\pm 20\%$

Emission Frequency

400KHz, 470KHz, 700KHz, 1000KHz $\pm 10\%$

Treatment modality

Capacitive and resistive

Capacitive heads

3 \varnothing sizes 30, 50, 70 mm

Resistive heads

3 \varnothing sizes mm 30, 50, 70 mm

Power supply: mains

100-240V AC, 50-60Hz, 1.8A

Size

170x220x60 mm


Weight

1100 gr



STANDARD EQUIPMENT

- 1 Professional Carrying case
- 1 TECARVET 4000 device
- 1 Counter electrode (plate)
- 1 Monopolar resistive handpiece with 3 heads, \varnothing mm 30-50-70
- 1 Monopolar capacitive handpiece with 3 heads, \varnothing mm 30-50-70
- 1 Cream flask 250 ml
- 3 Elastic bands
- 1 Power supply unit
- 1 Operating manual



- RESISTIVE
- CAPACITIVE
- BIPOLEAR
- MULTIUSER
- PROGRAMMABLE
- DISPLAY TOUCH SCREEN
- MAINS POWER SUPPLY

TECARVET 4000

TECAR THERAPY

PROGRAM LIST

DOG / CAT

Hematomas / Edemas

Osteoarthritis Capacitive

Osteoarthritis Resistive

Osteoarthritis Capacitive

Osteoarthritis Resistive

Acute muscle contracture

Chronic muscle contractures

Acute Dorsalgia / Lumbago

Chronic Dorsalgia / Lumbago

Spondylarthrosis

HORSE

Hematomas / Edemas

Osteoarthritis Capacitive

Osteoarthritis Resistive

Osteoarthritis Capacitive

Osteoarthritis Resistive

Acute muscle contracture

Chronic muscle contractures

Acute Dorsalgia / Lumbago

Chronic Dorsalgia / Lumbago

Spondylarthrosis

Tendon injuries

Trigger points