

TECAR BEAUTY 7000 MED

AESTHETIC TECARTHERAPY

THE BEST DEVICE WITH MEDICAL CERTIFICATION IN THE GLOBUS TECAR BEAUTY LINE. THE DEVICE IS EQUIPPED WITH A WIDE RANGE OF PROGRAMS INCLUDING BOTH AESTHETIC TREATMENTS AND PROTOCOLS FOR DIFFERENT MUSCULOSKELETAL PATHOLOGIES. THE TOUCHSCREEN NAVIGATION MENU ON A COLOR SCREEN INCLUDES AN ATLAS WITH SEVERAL PICTURES AND A PRACTICAL DESCRIPTION OF HOW TO PERFORM THE TREATMENTS. MOREOVER, TECAR BEAUTY 7000 MED ALLOWS EXPERT USERS TO CREATE AND STORE NEW PROGRAMS IN JUST A FEW STEPS.

TECHNICAL SPECIFICATIONS

Display

Size 7" color touch screen

Nominal Power (NP)

350W \pm 10%

Emission frequencies

300kHz, 400kHz, 448kHz, 470kHz, 700kHz,
1000kHz, 1200kHz, \pm 10%

Treatment methods

Capacitive and Resistive

Capacitive Electrodes

3 \varnothing (30, 50, 70 mm)

Resistive Electrodes

3 \varnothing (30, 50, 70 mm)

Mains power supply

- 230V AC, 50-60Hz, 3.7A
- 115V AC, 50-60Hz, 7.4A

Size

515x340x170 mm

Weight

11.5 Kg

STANDARD EQUIPMENT

· Cod. G6244

· 1 TECAR BEAUTY BEAUTY 7000 MED device

· 1 Counter electrode (flexible plate)

· 1 Monopolar Resistive handpiece with 3 electrodes
 \varnothing (30, 50, 70 mm)

· 1 Monopolar Capacitive handpiece with 3 electrodes
 \varnothing (30, 50, 70 mm)

· 1 Bipolar handpiece \varnothing (30 mm)

· 1 Hands-free treatment kit

· 1 1000 ml conductive cream bottle

· 1 Power supply cable

· 1 Operating manual



TECAR BEAUTY 7000 MED

AESTHETIC TECARTHERAPY

MEDICAL TREATMENT LIST

Rev. En_10_2024

| | | |
|--|---------------------------------|---|
| Hallux valgus | Plantar fasciitis | Disc protrusion |
| Arthritis | Phlebitis | Acute osteitis pubis |
| Arthrosis | Fracture acute phase | Chronic osteitis pubis |
| Cervical arthrosis | Fracture sub-acute phase | Rhizartriosis |
| Dorsal or lumbar arthrosis | Gonarthrosis | Subacromial impingement syndrome |
| Bursitis | 1st degree muscular injury | Piriformis syndrome |
| Adhesive capsulitis | 2nd degree muscular injury | Iliotibial band syndrome |
| Acute cervical pain | 3rd degree muscular injury | Ulnar styloiditis |
| Chronic cervical pain | 1st degree ligament injury | Acute Achilles tendinitis |
| Cervico-brachialgy | 2nd degree ligament injury | Chronic Achilles tendinitis |
| Baker's cyst | 3rd degree ligament injury | Acute tendinitis |
| Whiplash | Chronic lymphedema | Chronic tendinitis |
| Nerve compression | Post-traumatic lymphedema | Pes aserinus tendinitis |
| Patello-femoral chondropathy (cartilage lesion) | Acute Low Back Pain | De Quervain's tendinitis/Radial styloiditis |
| Muscle contracture | Chronic Lumbago | Patellar tendinitis |
| Coxarthrosis | Lumbosciatica | Tenosynovitis |
| Muscle cool down | Metatarsalgia | Sprain injuries with edema |
| Epicondylitis | Osteochondrosis/Osteochondritis | Sprain injuries without edema |
| Epitrochleitis | Periarthritis | Carpal tunnel |
| Disc herniation | Periostitis | |
| | Post-ACL surgery | |















• The **Pre-set Area** includes 26 default programs intended for specific pathologies and connected with 5 different treatment categories.

| TREATMENT CATEGORIES | SPECIFIC PATHOLOGIES | TREATMENT CATEGORIES | SPECIFIC PATHOLOGIES |
|---|----------------------|---------------------------------------|----------------------|
| Lymphatic and vascular pathologies | | Bone and articular pathologies | |
| • Chronic lymphedema | 1 | • Arthritis | 1 |
| • Post-traumatic lymphedema | 1 | • Arthrosis | 1 |
| Muscular and soft-tissue pathologies | | • Disc herniation | 1 |
| • Scars | 1 | • Fracture acute phase | 1 |
| • Muscle contracture | 1 | • Fracture sub-acute phase | 1 |
| • Muscle cool down | 1 | • Disc protrusion | 1 |
| • Fibrosis | 1 | • Sprain injuries | 1 |
| • 1st degree muscular injury | 1 | Tendon and ligament injuries | |
| • 2nd degree muscular injury | 1 | • Ligament injury 1st degree | 1 |
| • 3rd degree muscular injury | 1 | • Ligament injury 2nd degree | 1 |
| Nervous pathologies | | • Ligament injury 3rd degree | 1 |
| • Cervico-brachialgy | 1 | • Acute tendinitis | 1 |
| • Nerve compression | 1 | • Chronic tendinitis | 1 |
| • Lumbosciatica | 1 | | |
| • Piriformis syndrome | 1 | | |
| • Carpal tunnel | 1 | | |

TECAR BEAUTY 7000 MED

AESTHETIC TECARTHERAPY

ESTHETIC ATLAS PROGRAM LIST

| | | |
|---|--|--|
| <p>FACE</p> <hr/> <p>Wrinkles</p> <hr/> <p>Eye contour</p> <hr/> <p>Acne scars</p> <hr/> <p>Lifting</p> <hr/> <p>Rosacea</p> <hr/> <p>Relaxing massage</p> <p>CLEAVAGE</p> <hr/> <p>Face-lift</p> <hr/> <p>Scars</p> <hr/> <p>Firming</p> <p>ARM-HAND</p> <hr/> <p>Wrinkles</p> <hr/> <p>Firming</p> <hr/> <p>Rejuvenation</p> <hr/> <p>Scars</p> <hr/> <p>Stretch marks</p> <hr/> <p>Fat reduction</p> <p>ABDOMEN-HIPS</p> <hr/> <p>Fat reduction</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Stretch marks</p> | <p>GLUTEUS</p> <hr/> <p>Stretch marks</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Lymphatic drainage</p> <hr/> <p>Fat reduction</p> <hr/> <p>Edematous cellulite</p> <hr/> <p>Fibrous cellulite</p> <hr/> <p>Sclerotic cellulite</p> <p>Massage</p> <p>BACK</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Massage</p> <hr/> <p>Fat reduction</p> <p>THIGHS AND LEGS</p> <hr/> <p>Stretch marks</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Lymphatic drainage</p> <hr/> <p>Fat reduction</p> <hr/> <p>Edematous cellulite</p> <hr/> <p>Fibrous cellulite</p> <hr/> <p>Sclerotic cellulite</p> <hr/> <p>Massage</p> | <p> TREATMENT ATLAS</p> <p> PULSED TREATMENTS</p> <p> ATHERMAL TREATMENTS</p> <p> RESISTIVE</p> <p> CAPACITIVE</p> <p> MULTI-FREQUENCY</p> <p> AUTOMATIC SCANNING OF FREQUENCIES</p> <p> MYOFASCIALS</p> <p> TRANSDERMAL DELIVERY</p> <p> HANDS FREE MODE</p> <p> BIPOLAR</p> <p> COLOR DISPLAY</p> <p> TOUCH SCREEN DISPLAY</p> <p> MAINS POWER SUPPLY</p> |
|---|--|--|

The **Pre-set/Aesthetic Area** includes 13 default programs intended for specific pathologies and connected with 5 different treatment categories.

| TREATMENT CATEGORIES | PROGRAMS | TREATMENT CATEGORIES | PROGRAMS |
|----------------------------|----------|-----------------------------------|----------|
| Lymphatic drainage | 1 | Lassità cutanea | |
| Localized fat | | • Minor | 1 |
| • Fold < 4 cm | 1 | • Moderate | 1 |
| • Fold 4-8 cm | 1 | • Serious | 1 |
| • Fold > 8 cm | 1 | | |
| Cellulite treatment | | Face-cleavage rejuvenation | |
| • Minor cellulite | 1 | • Minor | 1 |
| • Moderate cellulite | 1 | • Moderate | 1 |
| • Serious cellulite | 1 | • Serious | 1 |

