

TECAR BEAUTY 7000

AESTHETIC TECARTHERAPY

POWERFUL AND VERSATILE DEVICE FROM THE TECAR BEAUTY RANGE. NAVIGATION SYSTEM WITH COLOR TOUCHSCREEN TO MAKE TREATMENT SEARCH AND PARAMETER MODIFICATION EASIER, AN ATLAS THAT DESCRIBES THE TREATMENTS DIVIDED INTO DIFFERENT PARTS OF THE BODY WITH COLOR PHOTOS AND A SECTION WHERE THE OPERATOR CAN STRUCTURE AND STORE PERSONALIZED TREATMENTS.
THE BEST EXPRESSION OF TECAR BEAUTY FOR BEAUTY TREATMENTS!

TECHNICAL SPECIFICATIONS

Display

Size 7" color touch screen

Power

350W \pm 10%

Emission frequencies

400kHz, 448kHz, 470kHz, 700kHz,
1000kHz, 1200kHz, \pm 10%

Treatment methods

Capacitive and Resistive

Capacitive Electrodes

3 \varnothing (30, 50, 70 mm)

Resistive Electrodes

3 \varnothing (30, 50, 70 mm)

Mains power supply

- 230V AC, 50-60Hz, 3.7A
- 115V AC, 50-60Hz, 7.4A

Size

515x340x170 mm

Weight

11.5 Kg

STANDARD EQUIPMENT

· Cod. G6248

· 1 TECAR BEAUTY 7000 device

· 1 Counter electrode (flexible plate)

· 1 Monopolar Resistive handpiece with 3 electrodes
 \varnothing (30, 50, 70 mm)

· 1 Monopolar Capacitive handpiece with 3 electrodes
 \varnothing (30, 50, 70 mm)

· 1 Bipolar handpiece \varnothing (30 mm)

· 1 Hands-free treatment kit

· 1 500 ml hyaluronic acid conductive face cream

· 1 1000 ml conductive cream for lymphatic-drainage
with eascin, fucus and caffeine

· 1 Power supply cable

· 1 Operating manual

CE

















GLOBUS
ITALIAN EXCELLENCE

TECAR BEAUTY 7000

AESTHETIC TECARTHERAPY

ESTHETIC ATLAS PROGRAM LIST

<p>FACE</p> <hr/> <p>Wrinkles</p> <hr/> <p>Eye contour</p> <hr/> <p>Acne scars</p> <hr/> <p>Lifting</p> <hr/> <p>Rosacea</p> <hr/> <p>Relaxing massage</p> <p>CLEAVAGE</p> <hr/> <p>Face-lift</p> <hr/> <p>Scars</p> <hr/> <p>Firming</p> <p>ARM-HAND</p> <hr/> <p>Wrinkles</p> <hr/> <p>Firming</p> <hr/> <p>Rejuvenation</p> <hr/> <p>Scars</p> <hr/> <p>Stretch marks</p> <hr/> <p>Fat reduction</p> <p>ABDOMEN-HIPS</p> <hr/> <p>Fat reduction</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Stretch marks</p>	<p>GLUTEUS</p> <hr/> <p>Stretch marks</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Lymphatic drainage</p> <hr/> <p>Fat reduction</p> <hr/> <p>Edematous cellulite</p> <hr/> <p>Fibrous cellulite</p> <hr/> <p>Sclerotic cellulite</p> <p>Massage</p> <p>BACK</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Massage</p> <hr/> <p>Fat reduction</p> <p>THIGHS AND LEGS</p> <hr/> <p>Stretch marks</p> <hr/> <p>Firming</p> <hr/> <p>Scars</p> <hr/> <p>Lymphatic drainage</p> <hr/> <p>Fat reduction</p> <hr/> <p>Edematous cellulite</p> <hr/> <p>Fibrous cellulite</p> <hr/> <p>Sclerotic cellulite</p> <hr/> <p>Massage</p>	<p> TREATMENT ATLAS</p> <p> PULSED TREATMENTS</p> <p> ATHERMAL TREATMENTS</p> <p> RESISTIVE</p> <p> CAPACITIVE</p> <p> MULTI-FREQUENCY</p> <p> AUTOMATIC SCANNING OF FREQUENCIES</p> <p> MYOFASCIALS</p> <p> TRANSDERMAL DELIVERY</p> <p> HANDS FREE MODE</p> <p> BIPOLAR</p> <p> COLOR DISPLAY</p> <p> TOUCH SCREEN DISPLAY</p> <p> MAINS POWER SUPPLY</p>
---	--	--

The **Pre-set/Aesthetic Area** includes 13 default programs intended for specific pathologies and connected with 5 different treatment categories.

TREATMENT CATEGORIES	PROGRAMS	TREATMENT CATEGORIES	PROGRAMS
Lymphatic drainage	1	Lassità cutanea	
Localized fat		• Minor	1
• Fold < 4 cm	1	• Moderate	1
• Fold 4-8 cm	1	• Serious	1
• Fold > 8 cm	1		
Cellulite treatment		Face-cleavage rejuvenation	
• Minor cellulite	1	• Minor	1
• Moderate cellulite	1	• Moderate	1
• Serious cellulite	1	• Serious	1

