

TECAR BEAUTY ELITE

AESTHETIC TECARTHERAPY

TECAR BEAUTY IS A NON-INVASIVE METHOD THAT FIGHTS AGING AND SKIN LAXITY. THE ELECTROMAGNETIC WAVES PRODUCED BY THE DEVICE PENETRATE INTO THE DERMIS AND RELEASE HEAT IN THE TISSUES. THIS CAUSES THE SYNTHESIS OF NEW COLLAGEN AND ELASTIN, VISIBLY IMPROVING THE COMPACTNESS OF THE SKIN. IN ADDITION, THE HEAT FACILITATES MICROCIRCULATION, PROMOTING A BETTER OXYGENATION OF THE TISSUES. PROFESSIONAL DEVICE SUITABLE FOR AESTHETIC CENTERS.

TECHNICAL SPECIFICATIONS

Display

Size 7" color touch screen

Power

50W ±10% Capacitive

25W ±10% Resistive

Emission frequencies

400kHz, 448kHz, 470kHz, 700kHz,

1000kHz, 1200kHz, ± 10%

Treatment methods

Capacitive and Resistive

Capacitive Electrodes

3 Ø (30, 50, 70 mm)

Resistive Electrodes

3 Ø (30, 50, 70 mm)

Mains power supply

- 230V AC, 50-60Hz, 3.7A

- 115V AC, 50-60Hz, 7.4A

Size

515x340x170 mm

Weight

11.5 Kg

STANDARD EQUIPMENT

· Cod. G6243

· 1 TECAR BEAUTY ELITE device

· 1 Counter electrode (flexible plate)

· 1 Monopolar Resistive handpiece with 3 electrodes
Ø (30, 50, 70 mm)

· 1 Monopolar Capacitive handpiece with 3 electrodes
Ø (30, 50, 70 mm)

· 1 Bipolar handpiece Ø (30 mm)

· 1 Hands-free treatment kit

· 1 500 ml hyaluronic acid conductive face cream

· 1 1000 ml conductive cream for lymphatic-drainage
with eascin, fucus and caffeine

· 1 Power supply cable

· 1 Operating manual



In compliance with Interministerial Decree no. 206/2015, the device is suitable for beauty centers.



TECAR BEAUTY ELITE

AESTHETIC TECARTHERAPY

ESTHETIC ATLAS PROGRAM LIST

FACE

Wrinkles

Eye contour

Acne scars

Lifting

Rosacea

Relaxing massage

CLEAVAGE

Face-lift

Scars

Firming

ARM-HAND

Wrinkles

Firming

Rejuvenation

Scars

Stretch marks

Fat reduction

ABDOMEN-HIPS

Fat reduction

Firming

Scars

Stretch marks

GLUTEUS

Stretch marks

Firming

Scars

Lymphatic drainage

Fat reduction

Edematous cellulite

Fibrous cellulite

Sclerotic cellulite

Massage

BACK

Firming

Scars

Massage

Fat reduction

THIGHS AND LEGS

Stretch marks

Firming

Scars

Lymphatic drainage

Fat reduction

Edematous cellulite

Fibrous cellulite

Sclerotic cellulite

Massage



TREATMENT ATLAS



RESISTIVE



CAPACITIVE



MULTI-FREQUENCY



MYOFASCIALS



HANDS FREE MODE



BIPOLAR



COLOR DISPLAY



TOUCH SCREEN DISPLAY



MAINS POWER SUPPLY

The **Pre-set/Aesthetic Area** includes 13 default programs intended for specific pathologies and connected with 5 different treatment categories.

TREATMENT CATEGORIES	PROGRAMS	TREATMENT CATEGORIES	PROGRAMS
Lymphatic drainage	1	Lassità cutanea	
Localized fat		• Minor	1
• Fold < 4 cm	1	• Moderate	1
• Fold 4-8 cm	1	• Serious	1
• Fold > 8 cm	1		
Cellulite treatment		Face-cleavage rejuvenation	
• Minor cellulite	1	• Minor	1
• Moderate cellulite	1	• Moderate	1
• Serious cellulite	1	• Serious	1



DOMINO S.R.L.
Via Vittorio Veneto, 52
31013 Codognè (TV) - ITALY



(+39) 0438 7933



globuscorporation.com

GLOBUS
ITALIAN EXCELLENCE