

MAGNUM 3000 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

POWERFUL PROFESSIONAL MAGNETOTHERAPY EQUIPPED WITH MANY SPECIAL PROGRAMS AND FUNCTIONS THAT MAKE IT EFFECTIVE THE TREATMENT OF SEVERAL OSTEOARTICULAR PATHOLOGIES AND THE RICOVERY AFTER ORTHOPEDIC SURGERY. POSSIBILITY OF LOCAL AND TOTAL-BODY TREATMENTS.

TECHNICAL SPECIFICATIONS

Display

Size of the visible area 2.6"

Channels

2 independent

Intensity can be emitted with flexible diffuser

- 200 Gauss on average for channel
- 400 Gauss on average on the whole (for 2 channels)

Frequency

5-200 Hz

Duty Cycle

5-50%

Free programming

Possibility of adding up to 10 custom programs

1+1 function

Possibility of running 2 different programs at the same time, using 1 channel each

Power supply

- Mains
- Rechargeable batteries (not included)


Size | Weight

170x220x60 mm | 820 gr

STANDARD EQUIPMENT

Code **G5335**

- 1 Bag
- 1 MAGNUM 3000 PRO magnetotherapy device
- 2 Flexible solenoids, 300x100 mm
- 1 Magnum Test
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTION

Code **G6034** MAGNUM 3000 PRO with 2 solenoids Pocket Pro



CE
0476



- EASY PROGRAM**
- 1+1** 1+1 MODE
- Last 10** LAST 10
- PRO PROGRAMMABLE**
- MULTI USER**
- RENT** RENTAL OPTION
- STOP RUN TIME** RUN TIME / EDIT PARAMETERS
- PLUG** MAINS SUPPLY

GLOBUS
ITALIAN EXCELLENCE

MAGNUM 3000 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

PROGRAM LIST

ARTHROSIS/ARTHRITIS

Knee arthrosis

Cervical arthrosis

Arthritis

Arthrosis

Hand arthrosis

Hip arthrosis

Shoulder arthrosis

Column arthrosis

PAIN

Algodystrophy - Sudek disease

Back pain

Joint pain

INFLAMMATIONS

Anti-inflammatory

Tendinitis

Epicondylitis

Epitrochleytis

Carpal tunnel

Shoulder tendinitis

ORTHOPEDICS

Post ACL surgery

Cartilage lesions

Knee posthesis

Hip prosthesis

Spinal fusion

Osteotomy

BONE PATHOLOGIES

Fractures

Scaphoide fracture

Wrist fracture

Shoulder fracture

Collarbone fracture

Forearm fracture

Pelvic fracture

Tibial fracture

Femoral fracture

Delayed fracture consolidation

Fracture with plaster

Osteonecrosis

Osteoporosis

Pseudoarthrosis

Bone edemas

TRAUMAS / EDEMAS / WOUNDS

Ankle sprain

Edemas- contusions

Ulcers

NIGHT PROGRAMS

Night program 15 Hz

Night program 50 Hz

Night program 75 Hz

DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY

Multi-frequency program

Program 5 Hz

Program 10 Hz

Program 15 Hz

Program 20 Hz

Program 25 Hz

Program 30 Hz

Program 35 Hz

Program 40 Hz

Program 45 Hz

Program 50 Hz

Program 55 Hz

Program 60 Hz

Program 65 Hz

Program 70 Hz

Program 75 Hz

Program 80 Hz

Program 85 Hz

Program 90 Hz

Program 95 Hz

Program 100 Hz

Program 120 Hz

Program 140 Hz

Program 160 Hz

Program 180 Hz

Program 200 Hz