

# MAGNETOVET 4000 HORSE

## MAGNETOTHERAPY FOR HORSES

THE NEW MAGNETOVET 4000 HORSE PRO KIT, IN ADDITION TO THE DEVICE FOR MAGNETOTHERAPY WITH THE NEW PROGRAMS DEDICATED TO THE HORSE, INCLUDES ACCESSORIES DESIGNED FOR ITS TREATMENT: THERAPEUTIC BLANKET AND NECK-COVER, FLEXIBLE SOLENOIDS WITH RELATED SPLITTER CABLES, A PRACTICAL DEVICE CASE AND AN ADJUSTABLE STRAP FOR ATTACHING THE DEVICE TO THE HORSE. A CONVENIENT CARRYING BAG AND A DOUBLE BATTERY PACK COMPLETE THE KIT. MAGNETOVET 4000 HORSE IS THE HIGH INTENSITY PULSED ELECTROMAGNETIC FIELDS MAGNETOTHERAPY DEVICE SPECIALLY DESIGNED FOR THE TREATMENT OF HORSES. THE FOUR CHANNELS ALLOW THE SOLENOIDS, WHICH ARE INSERTED INTO THE THERAPEUTIC BLANKET AND NECK-COVER, TO BE CONNECTED. THIS ALLOWS THE ENTIRE SPINE OF THE HORSE TO BE TREATED.

### TECHNICAL FEATURES

#### Display

Size 2.6", backlit

#### Channels

4

#### Intensities

200 Gauss on average per channel, with flexible solenoids

800 Gauss on average as a whole (on 4 channels)

The nominal powers refer to accessories with two coils. In the case of accessories with multiple coils, the power is evenly distributed

#### Frequency

Adjustable in programming from 5 to 200Hz

#### Duty Cycle

5% a 50%. Increment by steps of 5

#### Power supply

- Main

- Rechargeable batteries 4000mAh included

#### Size - Weight

17x22x6 cm - 1,37 Kg

### STANDARD EQUIPMENT


Cod. G6589

· 1 Bag

· 1 MAGNETOVET 4000 HORSE device

· 1 Power supply unit

· 1 Operating manual



CE



- 2+2 MODE
- LAST 10
- PROGRAMMABLE
- BATTERY PACK
- MAINS POWER SUPPLY

# MAGNETOVET 4000 HORSE

## MAGNETOTHERAPY FOR HORSES

Rev. En\_04\_2024

### LIST OF PROGRAMS FOR HORSE

THE PROGRAMS ARE DIVIDED BY THE ANATOMICAL AREAS OF THE NECK, BACK AND LIMBS

Arthrosis

---

Arthritis

---

Fracture

---

Delayed fracture consolidation

---

Discopathy

---

Tendinopathy

---

Wound

---

Ulcer

---

Anti-inflammatory

---

Injury

---

Contusion

---

Edema

---

Muscular relaxation

---

Decontracting

---

Muscular warm-up

---

Post-competition recovery

---

Muscle injury

### LIST OF PROGRAMS FOR DOG/CAT

THE PROGRAMS INVOLVE THE USE OF THE THERAPEUTIC BED OR FLEXIBLE SOLENOIDS WITH 2 COILS, 30x10 cm

Fracture

---

Osteoarthritis

---

Discopathy

---

Cicatrization

---

Inflammatory pain