

DERMATO-FUNCTIONAL TECAR

A single technology, **infinite application possibilities**



GLOBUS
ITALIAN EXCELLENCE



ABOUT US

GLOBUS has been a globally successful company producing portable electromedical equipment for over 30 years.

Our wide range includes complete lines of products for **PHYSIOTHERAPY, SPORTS, VETERINARY** and **AESTHETICS**.

OUR MISSION

The mission of **GLOBUS** is to contribute to promoting health and well-being in the world by providing high quality products for rehabilitation and physiotherapy professionals as well as to private users who wish to improve their quality of life in total safety and benefit from the most modern technologies available on the market.

OUR HISTORY

GLOBUS was founded in 1984 in Codognè, in the Province of Treviso (Italy), by its current President and CEO, Pierpaolo Lucchetta, bronze in volleyball at the Los Angeles Olympics and award-winning professional until 1994.

The high-level athlete experience of its President has been fundamental for **GLOBUS** to design and develop different product lines for post-injury recovery, muscle strengthening and the general improvement of sport performance. In order to develop safe, effective and highly technological products, **GLOBUS** collaborates with several universities and research institutes all over the world and scrupulously adheres to the most rigorous scientific guidelines.

Today, **GLOBUS** is a renowned company in the rehabilitation and physiotherapy sector and exports to more than 50 countries around the world.

GLOBUS products are certified UNI CEI EN ISO 13485:2016 and Quality Assurance System CE in Europe, FDA in the United States, and IEx in Brazil.

TECAR GLOBUS Line

With the devices of the **TECAR GLOBUS** line, you can perform any type of **dermato-functional** treatment and work on multiple districts at all levels in an **effective, safe** and **rapid way**.

The **treatment guide**, the **atlas**, the **pre-set protocols**, the **special functions** and the wide range of **accessories** specifically developed for different application areas guarantee new possibilities and **always satisfactory results**.

REACTIVATE

CELLULAR METABOLISM

TECAR THERAPY acts with extraordinary effectiveness in depth on cellular metabolism and allows to obtain incredible results in different pathologies, including chronic ones, like no other technology is able to do.

Thanks to the versatility of our devices, it will be possible to intervene, for example, on tissue reactivation, sclerotic cellulite or oedemas and deep hematomas with excellent results from the first treatment.

- **VENOUS AND LYMPHATIC CIRCULATION DYSFUNCTIONS**
- **EDEMATOUS, FIBROUS AND SCLEROTIC CELLULITE**
- **SKIN LAXITY AND TISSUE HYPOTONIA**
- **LOCALIZED ADIPOSITY**
- **REABSORPTION OF OEDEMAS AND HEMATOMAS**



VENOUS AND LYMPHATIC CIRCULATION DYSFUNCTION



With **TECAR THERAPY** it is possible to **open the lymphatic stations** and **drain the accumulation of interstitial fluids**, facilitating the **reabsorption of swelling** in the limbs and the rebalancing of the venous and lymphatic system.

EDEMATOUS, FIBROUS AND SCLEROTIC CELLULITE

Treatment in Resistive mode with the **special handpieces (FASCIA TOOLS KIT)** for the **resolution of fibrous nodules** and **remodelling of the connective tissue**.



STIMULATE

SKIN HEALTH


With the **TECAR GLOBUS** devices, you can also stimulate the production of collagen and elastin in the skin in an exceptionally precise and controlled way, stimulating the natural process of tissue regeneration and exponentially reducing the recovery and healing time, among others, of post-surgical scars and keloids.

- **DERMATO-FUNCTIONAL TREATMENTS
FACE ANTI-AGING**
- **STIMULATION OF COLLAGEN AND
ELASTIN PRODUCTION**
- **PRE-SURGICAL PREPARATION**
- **OEDEMA AND PAIN POST AESTHETIC
AND RECONSTRUCTIVE SURGERY**
- **KELOIDS AND HYPERTROPHIC SCARS**
- **ULCERS AND BURNS**






STIMULATION OF COLLAGEN AND ELASTIN PRODUCTION



SKIN REGENERATED WITH TECAR GLOBUS



KELOIDS AND HYPERTROPHIC SCARS



Keloids are skin fibrotic lesions resulting from an alteration in the normal healing process of a skin wound. The combined treatment with **TECAR THERAPY** and **FASCIA TOOLS** allows both to prevent its formation, and to act effectively on the fibrotic tissue, for both aesthetic and functional results.

RESTORE

INTIMATE COMFORT

The accuracy of our **TECAR** devices allows to act effectively and in full safety and comfort, also to solve problems strictly related to the genital area.

Specific accessories have been specially designed for both internal and external treatments, which improve the health of the pelvic floor, both in men and women, restore vaginal tone and elasticity and guarantee excellent results in case of scars from childbirth.

- IMPROVEMENT OF PELVIC FLOOR TONE
- CHILDBIRTH SCARS
- ABDOMINAL DIASTASIS
- ENDOPELVIC MYOFASCIAL SYNDROMES
- PERINEAL REHABILITATION FOR MEN AND WOMEN
- VAGINAL LAXITY

IMPROVEMENT OF PELVIC FLOOR TONE

The **INTIMA KIT** is designed to carry out endopelvic treatments in total safety.

The capacitive mode allows a deep regeneration of the mucous membranes in maximum comfort.

In combination with traditional handpieces and the **HANDS-FREE** kits, it provides a complete system for tone enhancement and for all pelvic floor dysfunction.



CHILD BIRTH SCARS



Pregnancy leaves a complete physical and emotional transformation on the mother's body and sometimes signs, like external or internal scars that can induce pain and severe functional dysfunctions, even after some time.

GLOBUS TECAR THERAPY includes functions and special electrodes for the rehabilitation of all these dysfunctions.

ABDOMINAL DIASTASIS



Application of the **"HANDS-FREE"** method for the treatment of abdominal diastasis.

The adhesive electrodes enable a deep energy transfer that stimulates the activation and synthesis of neo collagen in the diastasis area, at the same time the therapist can work on respiratory or abdominal muscle re-education.

PROFESSIONAL DIACARE LINE

DIACARE 7000

- Medical Protocols
- Nominal Power (NP): 350W \pm 10%
- Display: 7" color touch screen
- Emission frequencies: 300kHz - 400kHz - 448kHz - 470kHz
700kHz - 1000kHz - 1200kHz, \pm 10%
- Mains power supply: - 230V AC, 50-60Hz, 3,7A
- 115V AC, 50-60Hz, 7,4A



DIACARE 5000

- Medical Protocols
- Nominal Power (NP): 180W \pm 10%
- Display: 4,3" color touch screen
- Emission frequencies: 300kHz - 400kHz - 448kHz
470kHz - 700kHz - 1000kHz, \pm 10%
- Mains power supply: 100-240V AC, 50-60Hz, 1,8A



PROFESSIONAL TECAR BEAUTY MED LINE

TECAR BEAUTY 7000 MED

- Aesthetic and Medical Protocols
- Nominal Power (NP): 350W \pm 10%
- Display: 7" color touch screen
- Emission frequencies: 300kHz - 400kHz - 448kHz - 470kHz
700kHz - 1000kHz - 1200kHz, \pm 10%
- Mains power supply: - 230V AC, 50-60Hz, 3,7A
- 115V AC, 50-60Hz, 7,4A



TECAR BEAUTY 6000 MED

- Aesthetic and Medical Protocols
- Nominal Power (NP): 180W \pm 10%
- Display: 4,3" color touch screen
- Emission frequencies: 300kHz - 400kHz - 448kHz
470kHz - 700kHz - 1000kHz, \pm 10%
- Mains power supply: 100-240V AC, 50-60Hz, 1,8A



CONVEX CAPACITIVE AND RESISTIVE ELECTRODES KIT



CONCAVE CAPACITIVE ELECTRODES KIT



ACTIVE NEUTRAL AND BIPOLAR HANDPIECES



RESISTIVE PLATE COVER KIT



DYNAMIC HANDLE



INTIMA KIT



HANDS-FREE PLUS KIT

FOR DIACARE 7000 AND TECAR BEAUTY 7000 MED



FOR DIACARE 5000 AND TECAR BEAUTY 6000 MED



NEUTRAL PLATES - FASCIA TOOLS KIT - HANDPIECES AND ELECTRODES



DEVICE CART



GLOBUS

ITALIAN EXCELLENCE

To know in detail the technical characteristics of our devices

visit the website

globuscorporation.com

or scan the QR Code



DOMINO SRL

Via Vittorio Veneto, 52 - 31013 Codognè (TV) - ITALY



(+39) 0438 7933

Information material for operators in the sector.

The images and the technical features in this brochure are not binding and can be changed by the manufacturer at any time without prior notice.

For further information about our products and treatments, to watch the videos about the products and to find our dealers, visit our website.